

# Tennant BioModulator Pro - Modes

Mode Name / Frequency	Features – Types of things the program can help with
Assess (R) 60 Hz	Measurement – used to measure the BioTerminal points and other areas of the body, gives a ( <i>voltage reaction</i> ) reading. <b>This is used with the electrode on the back of the device, placed firmly on skin.</b>
Ten-8™ (R) 24-30 Hz	Pain & Injuries of Muscle, Tissue, Bone and Inflammation. - <b>Used with BioTransducer or any electrodes</b>
Infinity™ 29-121 Hz	Organs and Overall General Health - <b>*Use this when you don't know what program to use.</b>
Automatic (R) 29-121 Hz	Access (Measurement) + 1 Minute Therapy (Infinity Mode) This is helpful when you are measuring yourself and want to see if you can improve low areas with a quick 1 minute therapy of Infinity. <b>Use with electrode on back of device.</b>
SCR 147 (R) 147 Hz	Scars - This program is used to work on scars where energy blockage can be found, for the purpose of good energy flow through scar tissue. <b>Use with BioTransducer and soreness oil</b>
BioTerminal® Step (BT Step) 8-35 Hz	This program will do all the BioTerminal points for 30 seconds each, and repeats. <b>Use with footplates and hand grips or BioTransducer (Can be used to balance body and emotions)</b>
BioTerminal® Crown (BT-CR) (R) 15.18 Hz	Used to offer more voltage to this BioTerminal Point- Bring it back in balance.
BioTerminal® Head (BT-HD) (R) 13.53 Hz	Used to offer more voltage to this BioTerminal Point- Bring it back in balance.
BioTerminal® Neck (BT-NK) (R) 12.05 Hz	Used to offer more voltage to this BioTerminal Point- Bring it back in balance.
BioTerminal® Chest (BT-CH) (R) 10.74 Hz	Used to offer more voltage to this BioTerminal Point- Bring it back in balance.
BioTerminal® Abdomen(BT-Ab)(R) 10.14 Hz	Used to offer more voltage to this BioTerminal Point- Bring it back in balance.
BioTerminal® Pelvic (BT-PV) (R) 9.01 Hz	Used to offer more voltage to this BioTerminal Point- Bring it back in balance.
BioTerminal® Base (BT-BA) (R) 8.02 Hz	Used to offer more voltage to this BioTerminal Point- Bring it back in balance.
Visual Acuity (VAstep) 4-1555 Hz and Infinity	Used to help with eye therapies. This program is used with the BioTransducer placed at the side of temple inline with the eye. <i>See BioTransducer Instructions for specifics on eye treatments.</i> <b>*This is also a great program for all areas of the body</b>
Tennant Pain Management (TEN PM) 0.5 - 304 Hz	Pain & Injuries of Muscle, Tissue, Bone and Inflammation. - <b>Used with BioTransducer or any electrodes</b>
Solfeggio Musical Notes 174, 285, 396, 417, 528, 639, 741, 852, 963 Hz	This program has the frequencies of musical notes, resonates well with all areas of the body. (is used plugged into footplates and hand grips and with the Transducer) Overall General Health - <b>* This is one of the best programs to choose, when you don't know what program to use.</b>
Harmonics (Schumann Waves) 7.83-33.8 Hz	Grounding, like walking barefoot. Use for correcting polarity and when you have been around a lot of EMF, computers, Wi-Fi - <b>Used with BioTransducer for Polarity correction but can also be used for grounding with electrodes.</b>

Mode Name / Frequency	Features
Delta Brain Waves <span style="float: right;">0.5 - 4.0 Hz</span>	Deepest meditation and deep restorative sleep, Dreamless sleep. Healing and Rejuvenation. Delta waves suspend external awareness (located at back of head) - <b>Occipital Lobe - Use with the BioTransducer</b>
Theta Brain Waves <span style="float: right;">3 - 8 Hz</span>	Trance or hypnotic state more susceptible to hypnosis, daydream, twilight state, dream state, REM sleep. Subconscious, creativity and intuition. Improved memory, inspiration and spiritual connection, Improved physical healing, Release of beneficial hormones related to health and longevity. Reduction of mental fatigue, anxiety, and stress. Overall bliss and contentment (located above the ear) <b>Temporal lobe - Use with the BioTransducer</b>
Alpha Brain Waves <span style="float: right;">7 - 12 Hz</span>	'Frequency Bridge' between our conscious thinking (Beta) and subconscious (Theta) mind. known to help calm you down and promote feelings of deeper relaxation. Positive Thinking, Aid overall mental coordination, calmness, alertness, mind/body integration and learning, deep studying (located at top of the head) – <b>Parietal Lobe - BioTransducer</b>
Beta Brain Waves <span style="float: right;">12 - 31.2 Hz</span>	Normal Wakeful State, productivity, concentration, and alertness improved logic, reasoning, and critical thinking (located at front of head-) <b>Frontal Lobe - Use with the BioTransducer</b>
Gamma Brain Waves <span style="float: right;">32 - 43 Hz</span>	High Performance, perception is heightened, peak mental state. Improve memory and processing of information from different brain areas. Pass information rapidly and quietly. - <b>Cerebellum</b> - Lower back of the head- <i>Research suggests that people with learning difficulties or impaired mental processing may not produce as many gamma waves.</i>
Chakra Step (CH step) <span style="float: right;">126 - 272 Hz</span>	Cycles between all Chakra frequencies every 30 seconds, and repeats. - <b>Use with the BioTransducer</b>
Chakra Root (CH Rt) (R) <span style="float: right;">194 Hz</span>	Treat with BioTransducer from the back-Located at Base between the legs (Perineum)
Chakra Sacrum (CH Sac) (R) <span style="float: right;">210 Hz</span>	Treat with BioTransducer from the back-Located two inches below the navel and center of the pelvis (backside)
Chakra Solar Plexus (CH Solrp) (R) <span style="float: right;">126 Hz</span>	Treat with BioTransducer from the back-Located along your spine and above the navel, in the upper abdomen (backside)
Chakra Heart (CH Ht) (R) <span style="float: right;">136 Hz</span>	Treat with BioTransducer from the back-Located back center of chest (backside)
Chakra Throat (CH Thr) (R) <span style="float: right;">141 Hz</span>	Treat with BioTransducer from the back-Located back of Neck
Chakra 3rd Eye (CH 3rd) (R) <span style="float: right;">221 Hz</span>	Treat with BioTransducer from back of head opposite of Forehead
Chakra Crown (CH Crn) (R) <span style="float: right;">172 Hz</span>	Treat with BioTransducer –Top of Head
Chakra 8th (CH 8th)(R) <span style="float: right;">272 Hz</span>	Treat with BioTransducer from back -Located arms length above your head (above the crown)
Blue Relax(R) <span style="float: right;">82 Hz</span>	Similar to Assess but provides higher frequency output. Measures the body and gives a reaction reading.
Modulate <span style="float: right;">139 - 147 Hz</span>	Acute pain and inflammation.
Stimulate <span style="float: right;">28 - 111 Hz</span>	Scenar (FM/Var) Pain and inflammation- modulates both power output and frequency to reduce accommodation.
Deep Stim <span style="float: right;">30 - 121 Hz</span>	Scenar (SW/Var) Pain and inflammation- multiple pulses and modulates both power output and frequency.
Blue Stim <span style="float: right;">22, 163, 22 Hz</span>	Scenar - Pain and inflammation- this mode outputs multiple pulses and modulates frequency.
Acute <span style="float: right;">121 Hz</span>	Scenar Russian technology – New Injury

Mode Name / Frequency	Features
Power 1 <span style="float: right;">57 - 122 Hz</span>	Toning face (when using on the face keep the power level below 35 and only treat for 1 - 2 minutes on each area, on clean dry skin) - <b>Most people use the Face Electrode for toning the face.</b>
Power Plus <span style="float: right;">57 - 122 Hz</span>	Toning of muscles (has modulation) - <b>Use electrode pads for toning large muscle groups.</b>
Power Advanced <span style="float: right;">57 - 122 Hz</span>	Toning Body (higher intensity of pulses) <b>When toning the larger muscle groups like the abs the electrode pads work well.</b>
Power Intense <span style="float: right;">57 - 122 Hz</span>	Toning Abs Muscles, Arms, Thighs ( higher intensity of pulses) <b>Use electrode pads for toning large muscle groups.</b>
PG 2500 <span style="float: right;">2300 - 2500 Hz</span>	Highest Hertz- Outer energy field. <b>Use with the Transducer.</b>
HGH <span style="float: right;">1725, 645, 1342 Hz</span>	Human Growth Hormone - <b>Most people use this with the Transducer.</b>
Repetitive Stimulate (R-Stim) <span style="float: right;">2500 Hz</span>	Used by Physical Therapists
Fast T (R) <span style="float: right;">500 Hz</span>	Used by Physical Therapists- Fast Twitch Muscle Fibers
Slow T (R) <span style="float: right;">250 Hz</span>	Used by Physical Therapists- Slow Twitch Muscle Fibers
Lung (LU) (R) <span style="float: right;">824 Hz</span>	Use by Acupuncturist-to treat meridian and to target organ for treatment - <b>Used with BioTransducer or any electrodes</b>
Pericardium (PC) (R) <span style="float: right;">530 Hz</span>	Use by Acupuncturist-to treat meridian and to target organ for treatment - <b>Used with BioTransducer or any electrodes</b>
Heart (HT) (R) <span style="float: right;">497 Hz</span>	Use by Acupuncturist-to treat meridian and to target organ for treatment - <b>Used with BioTransducer or any electrodes</b>
Small Intestine (SI) (R) <span style="float: right;">791 Hz</span>	Use by Acupuncturist-to treat meridian and to target organ for treatment - <b>Used with BioTransducer or any electrodes</b>
Triple Burner (TE) (R) <span style="float: right;">732 Hz</span>	Chinese Medicine - Upper - Heart, Lungs, Pericardium, , Middle -Spleen, Stomach, Gallbladder and Lower -Liver, Kidneys, Urinary Bladder, Small Intestine, Large Intestine and Uterus. <b>Use BioTransducer at - TB04 is located 'on transverse crease of dorsum of wrist in depression lateral to tendon of extensor digitorum communis.</b>
Large Intestine (LI) (R) <span style="float: right;">553 Hz</span>	Use by Acupuncturist-to treat meridian and to target organ for treatment - <b>Used with BioTransducer or any electrodes</b>
Spleen (SP) (R) <span style="float: right;">702 Hz</span>	Use by Acupuncturist-to treat meridian and to target organ for treatment - <b>Used with BioTransducer or any electrodes</b>
Liver (LR) (R) <span style="float: right;">442 Hz</span>	Use by Acupuncturist-to treat meridian and to target organ for treatment - <b>Used with BioTransducer or any electrodes</b>
Kidney (KI) (R) <span style="float: right;">608 Hz</span>	Use by Acupuncturist-to treat meridian and to target organ for treatment - <b>Used with BioTransducer or any electrodes</b>
Bladder (BL) (R) <span style="float: right;">667 Hz</span>	Use by Acupuncturist-to treat meridian and to target organ for treatment - <b>Used with BioTransducer or any electrodes</b>
Gallbladder (GB) (R) <span style="float: right;">583 Hz</span>	Use by Acupuncturist-to treat meridian and to target organ for treatment - <b>Used with BioTransducer or any electrodes</b>
Stomach (ST) (R) <span style="float: right;">471 Hz</span>	Use by Acupuncturist-to treat meridian and to target organ for treatment - <b>Used with BioTransducer or any electrodes</b>
Advanced Mode (AVA) <span style="float: right;">0.5 - 1565 Hz</span>	Where you can program additional frequencies to create additional therapy modes.
Memorized Mode (VAZ 4, 3, 2, 1)	Where you save the personalized modes, you created with the AVA program.
Emergency (ER) <span style="float: right;">24 - 30 Hz</span>	Use in an Emergency for snake and spider bites as taught in the conference. Double select middle button takes you instantly to maximum power. - <b>This is used with the electrode on the back of the device, placed firmly on skin.</b>